Instructions for Use

For Use Under Emergency Use Authorization (EUA) Only
For in vitro diagnostic use

CLIA COMPLEXITY: High (U.S.)

NAME AND INTENDED USE

The OraQuick™ Ebola 2.0 Rapid Antigen Test is an in vitro diagnostic single-use immunoassay intended for the presumptive qualitative detection of antigens from viruses within the *Orthoebolavirus* genus (such as *Orthoebolavirus zairense*, *Orthoebolavirus bundibugyoense*, *Orthoebolavirus sudanense*, and *Orthoebolavirus tairnense*); however, it does not differentiate among these viruses.

Testing with the OraQuick™ Ebola 2.0 Rapid Antigen Test must only be performed when public health authorities have determined the need for this test. Testing is limited to laboratories certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C. §263a, that meet the requirements to perform high complexity tests, or in similarly qualified non-U.S. laboratories and facilities adequately equipped, trained, and capable of testing for *Orthoebolavirus* infection. The test is intended for use by experienced personnel who have received specific training on the use of the OraQuick™ Ebola 2.0 Rapid Antigen Test, training in the correct use of recommended personal protective equipment (PPE), and expertise in infectious disease diagnostic testing, including the safe handling of clinical specimens potentially containing *Orthoebolavirus*. Testing for Ebola disease must be performed in accordance with current guidelines provided by the appropriate public health authorities that address appropriate biosafety conditions, interpretation of test results, and coordination of testing, results, and patient management.

The OraQuick™ Ebola 2.0 Rapid Antigen Test is intended for use in:

- Venipuncture whole blood and fingerstick whole blood specimens as an aid in the diagnosis of Ebola disease in patients suspected of and with signs and symptoms consistent with Ebola disease and who have epidemiological risk factors for *Orthoebolavirus* exposure.
- Cadaveric oral fluid collected from recently deceased individuals with epidemiological risk factors who are suspected to have had Ebola disease at the time of death. The test is intended as an aid in determining Ebola disease at the time of death to inform decisions regarding the safe handling of cadavers to prevent disease transmission.

The OraQuick™ Ebola 2.0 Rapid Antigen Test results are presumptive. Definitive identification of Ebola disease requires additional testing and confirmation procedures in consultation with public health and/or other appropriate authorities to whom reporting is required. The diagnosis of Ebola disease must be made based on history, signs, symptoms, exposure likelihood, and other laboratory evidence in addition to the identification of *Orthoebolavirus*. Positive results indicate the presence of *Orthoebolavirus* antigens; clinical correlation with patient history and other diagnostic information is necessary to determine infection status. Positive results do not rule out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease. Negative results do not preclude *Orthoebolavirus* infection and should not be used as the sole basis for treatment of patient management decisions. Negative results should be considered in the context of history, signs, symptoms, and exposure likelihood.

The OraQuick™ Ebola 2.0 Rapid Antigen Test is not intended for general *Orthoebolavirus* infection screening, such as airport screening or contract tracing of individuals at risk of exposure without signs or symptoms of infection. The OraQuick™ Ebola 2.0 Rapid Antigen Test is not intended for use as a diagnostic test for oral fluid samples of living individuals.

The OraQuick™ Ebola 2.0 Rapid Antigen Test is authorized for use under the U.S. Food and Drug Administration's Emergency Authorization.

SUMMARY AND EXPLANATION OF THE TEST

Ebola disease is a severe, often-fatal disease in humans and nonhuman primates that has appeared sporadically since its initial recognition in 1976. *Orthoebolavirus* species are in the *Filoviridae* family. There are four species of *Orthoebolavirus* affecting humans: *Orthoebolavirus zairense*, *Orthoebolavirus bundibugyoense*, *Orthoebolavirus sudanense*, and *Orthoebolavirus taiense*. The presence of *Orthoebolavirus* antigens indicates that the individual may be currently infected and capable of transmitting the virus.

The OraQuick™ Ebola 2.0 Rapid Antigen Test utilizes a sandwich capture lateral flow immunoassay method to detect *Orthoebolavirus* antigens. Ebola antigens are captured and visualized by colloidal gold labeled with Ebola antibodies generating a visible line in the test zone for a positive sample.

PRINCIPLES OF THE TEST

The OraQuick™ Ebola 2.0 Rapid Antigen Test is a manually performed, visually read immunoassay for the qualitative detection of *Orthoebolavirus* antigens in human venipuncture whole blood, fingerstick whole blood, and cadaveric oral fluid. The OraQuick™ Ebola 2.0 Rapid Antigen Test is comprised of both a single-use test device and vial containing a pre-measured amount of a buffered developer solution. The test consists of a sealed pouch with two separate compartments for each component. The OraQuick™ Ebola 2.0 Rapid Antigen Test utilizes a proprietary lateral flow immunoassay procedure.


The assay test strip, which can be viewed through the test device result window, is comprised of a series of components: the blocker pad, the conjugate pad, the nitrocellulose membrane, and finally the absorbent pad. The performance of the assay occurs through hydration and transport of reagents as they interact with the specimen across the strip via chromatographic lateral flow. As the specimen migrates through the device, *Orthoebolavirus* antigens, if present, are captured by antibodies located in the blocker pad and subsequently form complexes with antibody-labeled gold colorimetric reagent on the assay strip. These labeled complexes are immobilized at the test (T) zone, producing a visible purple line. In the absence of *Orthoebolavirus* antigens, no complexes form at the test (T) zone, and no line is observed. The intensity of the test line is not directly proportional to the concentration of the virus in the specimen.

Fifty microliters (50 µL) of a fingerstick or venipuncture whole blood specimen are transferred to the integrated collection pad on the device, followed by the insertion of the device into the developer vial. For cadaveric oral fluid testing, swab the oral cavity with the integrated with the integrated flat collection pad of the device and then insert the device directly into the developer solution to run the test. The developer solution facilitates the capillary flow of the specimen into the device and onto the assay strip. As the specimen flows through the device, antigens from the specimen are bound by the *Orthoebolavirus* antibody-labeled gold colorimetric reagent present on the assay strip. If the specimen contains *Orthoebolavirus* antigens, the resulting labeled complexes bind to the test (T) zone resulting in a purple line. If the specimen does not contain *Orthoebolavirus* antigens, the labeled complexes do not bind at the test zone and no line is observed. The intensity of the line color is not directly proportional to the amount of virus present in the specimen. The remaining colloidal gold is transported and bound to the control (C) zone. This procedural control serves to demonstrate that the fluid has migrated adequately through the device. A purple line will appear at the control (C) zone during the performance of all negative tests (refer to the *Test Result and Interpretation of Test Result* section on page 6 of in this package insert). Positive results may be interpreted as soon as purple lines are visible at the test (T) zone and control (C) zone. In instances of very high virus concentrations, the control line may not appear. If there is a test line present, the test should be interpreted as positive even if the control line is not present. Negative results MUST be read 30 to 40 minutes after inserting the device into the developer vial.

MATERIALS PROVIDED

The OraQuick™ Ebola 2.0 Rapid Antigen Test Kit is available in the following packaging configuration:

Components of Kit Catalog Number	25 count/kit 1001-0708_EUA
Divided pouch	25
Each pouch containing: Test device (1) Absorbent packet (1) Developer solution vial (1) (each vial contains 1.0 mL of a buffered saline solution with an antimicrobial agent)	
Test stands	25
Plastic capillary blood collection tubes	30
Instructions for use	1
Cadaveric oral fluid quick reference guide	1
Whole blood quick reference guide	1
Visual Reference card	1
Fact sheet for relatives and caregivers	1 (pack of 5)
Fact sheet for patients	1 (pack of 5)
Fact sheet for healthcare providers	1
Fact sheet for response team	1

MATERIALS REQUIRED AND AVAILABLE AS AN ACCESSORY TO THE KIT

OraQuick™ Ebola 2.0 Rapid Antigen Kit Controls 1001-0709_EUA

Ebola 2.0 Positive Control (1 vial, 0.5 mL); Ebola 2.0 Negative Control (1 vial, 0.5 mL); Package Insert

MATERIALS REQUIRED BUT NOT PROVIDED

Timer or watch capable of timing 30 to 40 minutes

Biohazard waste container

Materials required for venipuncture whole blood specimen collection or

Materials required to obtain a fingerstick whole blood specimen

Pipette capable of delivering 50 µL

Disposable pipette tips

WARNINGS AND PRECAUTIONS

- For prescription use only.
- This product has not been U.S. FDA cleared or approved, but has been authorized for emergency use by U.S. FDA under an EUA for use by the authorized laboratory.
- This product has been authorized only for the detection of antigens from orthoebolaviruses including *Orthoebolavirus zairense*, *Orthoebolavirus bundibugyoense*, *Orthoebolavirus sudanense*, and *Orthoebolavirus taiense*, not for any other viruses or pathogens.
- The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for the detection of Ebola virus under section 564(b) (1) of the Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.
- All testing **MUST** be conducted under appropriate biosafety conditions in accordance with applicable country, state and local laws and with CDC and/or WHO guidelines.
- In the U.S., Ebola disease is a nationally notifiable condition and must be reported to public health authorities in accordance with local, state, and federal regulations.
- Specimens **MUST** always be treated as infectious and/or biohazardous. The use of all possible universal precautions is highly recommended when handling specimens with this test.
- Use personal protective equipment (PPE) consistent with current guidelines including safety goggles and / or face shields, masks or respiratory equipment, disposable gowning, boots, and gloves. Users performing this test **MUST** be appropriately trained on the donning and doffing of personal protective equipment.
- All personnel conducting testing **MUST** read and be familiar with Universal Precautions¹, Infection Control for Viral Hemorrhagic Fevers in the African Health Care Setting² and in Information for Healthcare Worker in the United States (<https://www.cdc.gov/vhf/Ebola/healthcare-us/index.html>) depending upon their location of testing.
- All equipment and biohazardous waste **MUST** be discarded in accordance with country, state, and local laws and policies.
- This test kit is for use with venipuncture and fingerstick whole blood specimens and cadaveric oral fluid swab specimens only.
- Cadaveric oral swab specimens should not be stored and should be tested immediately after collection.
- If samples are suspected of hemolysis, another venous whole blood should be collected and tested.
- Do not freeze whole venous whole blood samples.
- Results must be read no earlier than 30 minutes and no later than 40 minutes.
- Read times below 30 minutes significantly affect performance of the device with low positive samples.
- Reducing the read time below 30 minutes can prevent correct identification of low positive samples.
- Reading any test result after 40 minutes may yield inaccurate test results.
- Samples **MUST NOT** be pre-diluted. Such pre-dilution can result in false negative results.
- This test is not intended to be used as a screening test on patients without signs and symptoms of Ebola infection or to monitor individuals who are undergoing treatment.
- Do not smoke, eat, or drink in areas where specimens or kit reagents are handled.
- This package insert must be read completely before using the product.
- Follow the instructions carefully when performing the OraQuick™ Ebola 2.0 Rapid Antigen Test. Failure to do so may cause an inaccurate test result.
- This test should be performed at temperatures in the range of 15°-40°C (59°-104°F). If stored refrigerated, ensure that the divided pouch is brought to operating temperature of 15°-40°C (59°-104°F) before performing testing.
- Do not use this test beyond the expiration date printed on the outer packaging. Always check expiration date prior to testing.
- Individuals who have received an Ebola vaccine may produce false results.

Device Handling Precautions

- Use all capillary tubes, pipette tips, test stands, test devices, and developer solution vials only once and dispose of properly (see *Warnings and Precautions*). **Do not reuse any test components.**
- Inspect the divided pouch. If the divided pouch has been damaged, discard the divided pouch and its contents and select a new divided pouch for testing.
- Do not interchange test devices and developer solution vials from kits with different lot numbers.
- Avoid microbial contamination and exercise care in handling the kit components.
- Adequate lighting is required to read a test result.

STORAGE INSTRUCTIONS AND SHELF LIFE

Store unused OraQuick™ Ebola 2.0 Rapid Antigen Tests unopened at 2°-30°C (36°-86°F). Do not open the divided pouch until you are ready to perform a test. If stored refrigerated, ensure that the divided pouch is brought to operating temperature 15°-40°C (59°-104°F) before opening.

DIRECTIONS FOR USE

GENERAL TEST PREPARATION

1. Follow the *Warnings and Precautions* section in this instructions for use.
2. Gather the materials you will need.
3. Allow the OraQuick™ Ebola 2.0 Rapid Antigen Test to come to operating temperature 15°-40°C (59°-104°F) before use.

- Set a test stand at your workspace, using only the stand provided.
- Open the two chambers of the divided pouch ("pouch") by tearing at the notches on the top of each side of the pouch (Picture 1, 2).
- Remove the developer solution vial ("vial") from the pouch. Hold the vial firmly in your hand. Carefully remove the cap from the vial by gently rocking the cap back and forth while pulling it off. Set the cap on your workspace cover.
- Slide the vial into the top of one of the slots in the Stand. **DO NOT** force the vial into the stand from the front of the slot as splashing may occur. Make sure the vial is pushed all the way to the bottom of the slot in the stand (Picture 3).
- Remove the device from the pouch. Handle the device by the plastic housing, **NOT** the collection pad (Picture 4). Place the device on a flat clean surface. Check to make sure that an absorbent packet is included with the device (Picture 5). If no absorbent packet is present, discard the device and obtain a new pouch for testing.

NOTE: DO NOT seal the two holes in the back of the device with labels or other materials. Doing so may cause an invalid result (Picture 6).

SPECIMEN COLLECTION AND TESTING PROCEDURE

The OraQuick™ Ebola 2.0 Rapid Antigen Test can be used for testing venipuncture whole blood and fingerstick whole blood specimens or cadaveric oral fluid specimens. Refer to the specific testing procedures below.

FINGERSTICK WHOLE BLOOD AND VENIPUNCTURE PROCEDURE

STEP 1: COLLECT THE SPECIMEN

STEP 1A: FINGERSTICK WHOLE BLOOD

- Using an antiseptic wipe, clean the finger of the person being tested. Allow the finger to dry thoroughly or wipe dry with a sterile gauze pad. Using a sterile lancet, puncture the skin just off the center of the finger pad and perpendicular to the fingerprint ridges. Apply slight pressure beside the point of the puncture (Picture 7). Avoid aggressive squeezing of the finger to make it bleed. Wipe away this first drop of blood with a sterile gauze pad. Allow a new drop of blood to form.
- Hold the capillary blood collection tube (Picture 8) horizontally and touch the tip of the tube to blood drop on the finger (Picture 9). The capillary tube fills automatically; do not squeeze the bulb during blood collection. To avoid air bubbles, keep the tube tip in continuous contact with the blood until the fill is complete. Confirm that the capillary tube is filled to the indicator line and that no bubbles are present (Picture 10). If bubbles are observed, discard the sample and collect a new one using a new capillary tube.

NOTE: If the capillary tube is dropped or comes in contact with any other surface, discard it in a biohazard waste container and obtain a new one. If blood is still flowing, attempt to collect 50 µL. If not, a new puncture may need to be performed. Get a new capillary tube for the collection of the blood sample.

STEP 1B: VENIPUNCTURE WHOLE BLOOD



- Using standard venous phlebotomy procedures collect a whole blood sample using a tube containing EDTA (lavender top) anticoagulant. If the specimens are not tested at the time of collection, the specimen may be stored at 2°-8°C (36°-46°F) for up to 6 days.
- Prior to testing, mix the blood tube gently by inversion several times to ensure a homogeneous sample.
- Using a calibrated pipettor, pipette 50 µL of sample and inspect the tip for bubbles (Picture 12). If a bubble is present, discard the collected sample and obtain a new sample.

STEP 2: TEST



- Deposit the blood sample onto the integrated collection pad at the position shown in the photos. (Pictures 13-14). **Samples MUST only be applied directly to the device and should not be pipetted into the developer vial.**

Note: Pipetting the sample at a location other than that shown in the photo could result in failure to detect a low positive sample.

Two chambers

DO NOT seal these two holes

2. Insert the collection pad of the device all the way into the developer vial (Picture 15). Make sure that the collection pad touches the bottom of the vial. The result window on the device should be facing you (Picture 16).
3. Start timing the test (Picture 17). **DO NOT** remove the device from the vial while the test is running. Fluid will appear and travel up the result window. The fluid may gradually lighten as the test develops (Picture 18). Read the results in a fully lit area. Read results between 30 and 40 minutes to obtain an accurate result. **DO NOT** read before 30 minutes or after 40 minutes.
4. Refer to the *Test Result and Interpretation of Test Result* section on page 6 of in this package insert.

CADAVERIC ORAL FLUID PROCEDURE – DIRECT COLLECTION

STEP 1: COLLECT THE SPECIMEN

1. Remove the device from its pouch. Handle the device by the plastic housing, **NOT** the collection pad (Picture 19). Check to make sure that an absorbent packet is included with the device (Picture 20). If no absorbent packet is present, discard the device and obtain a new pouch for testing.
2. Collect a sample from oral cavity by swiping the integrated flat collection pad across the upper and lower gum lines once (Pictures 21-22). Samples can also be collected by swabbing the cheek (Picture 23) or the soft palate (Picture 24). Only one sample location should be used when swabbing with the device.
NOTE: Both sides of the flat pad may be used during this procedure.

STEP 2: TEST

1. Insert the collection pad of the device all the way into the vial (Picture 25). Make sure that the collection pad touches the bottom of the vial. The result window on the device should be facing you (Picture 26).
2. Start timing the test (Picture 27). **DO NOT** remove the device from the vial while the test is running. Purple fluid will appear and travel up the result window. The purple fluid will gradually disappear as the test develops (Picture 28). Read results between 30 and 40 minutes to obtain an accurate result. **DO NOT** read before 30 minutes or after 40 minutes.
3. Refer to the *Test Result and Interpretation of Test Result* section on page 6 of this package insert.

NOTE: DO NOT cover the two holes in the back of the device with labels or other materials. Doing so may cause an invalid result (Picture 29).


GENERAL TEST CLEAN-UP

1. Dispose of the used test materials in a biohazard waste container. All equipment and biohazardous waste should be discarded in accordance with country, state, and local laws and policies.
2. Change your gloves between each test to prevent contamination. Throw away the used gloves in a biohazard waste container.
3. Use a freshly prepared 10% solution of bleach to clean up any spills.

QUALITY CONTROL PROCEDURES

Built-in Control Features

The OraQuick™ Ebola 2.0 Rapid Antigen Test has a built-in procedural control that demonstrates assay validity. A purple line in the control (C) zone of the result window indicates that the fluid migrated appropriately through the test device. In some patients with high levels of virus, a line in the C zone will be absent, but the line in the T zone will be positive. In this case, the T zone line serves as the built in procedural control. Refer to the *Test Result and Interpretation of Test Result* section on page 6 of this package insert.

Do NOT seal these two holes

External Quality Control

OraQuick™ Ebola 2.0 Rapid Antigen Test Kit Controls must be used with the OraQuick™ Ebola 2.0 Rapid Antigen Test. The kit controls are specifically formulated and manufactured to ensure performance of the Test, and are used to verify your ability to properly perform the test and interpret the results. The Ebola 2.0 positive control will produce a positive test result and the Ebola 2.0 negative control will produce a negative test result. Refer to the *Test Result and Interpretation of Test Result* section on page 6 of this package insert. Use of kit control reagents manufactured by any other source may not produce the required results, and therefore, may not meet the requirements for an adequate quality assurance program for the OraQuick™ Ebola 2.0 Rapid Antigen Test. If external controls do not produce expected results, patient testing should not be performed. Contact OraSure Technologies' Customer Care if the kit control reagents do not produce the expected result.

Run the External Controls under the following circumstances:

- Each new operator prior to performing testing on patient specimens,
- When opening a new test kit lot,
- When a new shipment of test kits is received,
- When test kit storage conditions fall outside of 2°-30°C (36°-86°F), and
- At periodic intervals as dictated by the user facility, country, state or local regulations and policies.

Test Procedure for External Controls:

Refer to the OraQuick™ Ebola 2.0 Rapid Antigen Test Kit Control package insert for full instruction on the use of these controls. It is the responsibility of each laboratory using the OraQuick™ Ebola 2.0 Rapid Antigen Test to establish an adequate quality assurance program to ensure the performance of the device under its specific locations and conditions of use.

TEST RESULT AND INTERPRETATION OF TEST RESULT

Interpret results 30 to 40 minutes after inserting the device into the developer vial. **Please refer to the Visual Reference Card provided in the kit to see colored images of the varying intensity levels that may be observed while running and interpreting the OraQuick™ Ebola 2.0 Rapid Antigen Test.**

NEGATIVE

A test is **Negative** if:

- A purple line appears in the C zone and NO line appears next to the T zone (Picture 30a/b)
- A negative test result is interpreted as Ebola antigen not detected in the specimen.
- The patient/cadaver is presumed negative for Ebola antigen.

A negative result does not preclude Ebola disease.

POSITIVE

A test is **Positive** if:

- A purple line appears in the T zone. Lines may vary in intensity. The test is positive regardless of how faint these lines appear (Picture 31a/b, 32a/b, 33a/b). In some cases the purple line in the C zone may not be present or may be very faint if there are high levels of virus in the sample (Picture 34a/b).
- A positive test result is interpreted as Ebola antigen detected in the specimen. The patient/cadaver is presumed positive for Ebola antigen. Definitive identification of Ebola disease requires additional testing and confirmation procedures in consultation with public health and/or other authorities to whom reporting is required.

In accordance with CDC and WHO recommendations, cadavers with a positive result should be subjected to safe and dignified burial procedures and contacts of an Ebola positive cadaver should be identified. Follow up testing of possible contacts should be conducted in accordance with, *EMERGENCY GUIDELINE Implementation and implementation of contact tracing for Ebola virus disease* (<https://iris.who.int/server/api/core/bitstreams/1be2b0a4-ecaa-442d-af2b-709dae0dae58/content>) issued by the World Health Organization (WHO) and the Centers for Disease Control (CDC).


INVALID

A test is **Invalid** if any of the following occurs:

- **NO** purple lines are present. (Picture 35), or
- A purple background in the result window makes it difficult to read the result at a time between 30 and 40 minutes (Picture 36), or
- Any partial line on one side of the C or T zones (Picture 37a/b, 38a/b).

An invalid test result means that there was a problem running the test, related either to the specimen or to the test device. **An invalid result cannot be interpreted. An invalid test result needs to be repeated with a fresh sample and a new device. Please contact OraSure Technologies' Customer Care if you are unable to obtain a valid test result upon repeat testing.**

To report a device problem, contact OraSure Technologies Customer Service (1-800-ORASURE) within the United States or +(001) 610-882-1820 for outside the United States. You can also go to <https://orasure.com/inquiry.html>.

LIMITATIONS OF THE TEST

1. U.S. hospitals or clinical laboratories concerned about a patient with potential *Orthoebolavirus* exposure should contact their local and/or state public health authorities. These agencies will work with CDC to determine whether a patient is or is not a person under investigation (PUI); only PUI should be tested for Ebola.
2. Testing with the OraQuick™ Ebola 2.0 Rapid Antigen Test should not be performed unless the patient meets clinical and epidemiologic criteria for testing suspect specimens.
3. All results with the OraQuick™ Ebola 2.0 Rapid Antigen Test are presumptive; definitive diagnosis of Ebola disease requires performing additional testing with molecular testing in consultation with public health and/or other authorities to whom reporting is required. The diagnosis of Ebola disease must be made based on history, signs, symptoms, exposure likelihood, and other laboratory evidence in addition to the identification of *Orthoebolavirus*.
4. Negative results do not preclude *Orthoebolavirus* infection and should not be used as the sole basis for patient management decisions.
5. Potential cross-reactivity of the OraQuick™ Ebola 2.0 Rapid Antigen Test with Ebola vaccines or therapeutics has not been evaluated. Specimens from patients who have received therapeutics or vaccines against *Orthoebolavirus* may exhibit false positive or other confounding test results.
6. The OraQuick™ Ebola 2.0 Rapid Antigen Test is a qualitative test; therefore, the intensity of the line for positive results should not be correlated with viral concentration.
7. The performance of this test was evaluated with a limited number of clinical specimens for whole blood, fingerstick blood, and cadaveric fluid.
8. The OraQuick™ Ebola 2.0 Rapid Antigen Test is not intended for use as a diagnostic test for oral fluid samples of living individuals nor is it intended for use with cadaveric samples stored in viral transport media.
9. False negative results have been observed with hemolyzed venous whole blood samples.
10. Testing of whole blood samples with concentrations of Rheumatoid Factor above 544 IU/mL may result in false positive results.
11. Cross-reactivity with organisms other than those tested in the cross-reactivity study have not been assessed.
12. The level of *Orthoebolavirus* antigens present in clinical specimens from individuals with early systemic infection is unknown. Test performance is associated with the level of antigen in the patient; therefore, the test is not intended for use in asymptomatic populations for mass screening purposes (e.g., as the sole means of Ebola disease control at airports or border crossings) or for testing individuals at risk of exposure without observable signs or symptoms of infection (e.g., contact tracing).

CONDITIONS OF AUTHORIZATION FOR THE LABORATORY

The OraQuick™ Ebola 2.0 Rapid Antigen Test, along with the authorized Fact Sheet for Healthcare Providers, the authorized Fact Sheet for Patients, the authorized Fact Sheet for Response Teams, the authorized Fact Sheet for Relatives and Caregivers and authorized labeling are available on the FDA website: <https://www.fda.gov/medical-devices/emergency-situations-medical-devices/emergency-use-authorizations-medical-devices#ebola>.

To assist clinical laboratories running OraQuick™ Ebola 2.0 Rapid Antigen Test, the relevant Conditions of Authorization are listed verbatim below, and are required to be met by laboratories performing the EUA test.

- Authorized laboratories* that receive your product must notify the relevant public health authorities of their intent to run your product prior to initiating testing.
- Authorized laboratories using your product must have a process in place for reporting test results to healthcare providers and relevant public health authorities, as appropriate.
- Authorized laboratories using your product must include with test result reports, all authorized fact sheets. Under exigent circumstances, other appropriate methods for disseminating these fact sheets may be used, which may include mass media.
- Authorized laboratories using your product must use your product as outlined in the authorized labeling. Deviations from the authorized procedures, including the authorized clinical specimen types, authorized control materials, authorized other ancillary reagents and authorized materials required to use your product are not permitted.
- Authorized laboratories must have a process in place to track adverse events and report to OraSure Technologies and to U.S. FDA pursuant to 21 CFR Part 803. For customers within the U.S., phone: 1-800-ORASURE (800-672-7873). For customers outside the U.S., phone +(001) 610-882-1820 or go to www.orasure.com.
- All laboratory personnel using your product must be appropriately trained in performing and interpreting the results of your product, use appropriate laboratory and personal protective equipment when handling your product, and use your product in accordance with the authorized labeling.
- OraSure Technologies, Inc., authorized distributor(s), and authorized laboratories must collect information on the performance of your product and must report any significant deviations from the established performance characteristics of your product of which they become aware to DMD/OHT/OPEQ/CDRH (via email: CDRH-EUA-Reporting@fda.hhs.gov). In addition, authorized distributor(s) and authorized laboratories report to OraSure Technologies. For customers within the U.S. phone: 1-800-ORASURE (800-672-7873). For customers outside the U.S., phone +(001) 610-882-1820 or go to www.orasure.com.
- OraSure Technologies, Inc., authorized distributor(s), and authorized laboratories using this product must ensure that any records associated with this EUA, are maintained until otherwise notified by U.S. FDA. Such records must be made available to U.S. FDA for inspection upon request.

*Authorized laboratories are laboratories certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C. §263a, and meet requirements to perform high-complexity tests, or in similarly qualified non-U.S. laboratories and facilities adequately equipped, trained, and capable of testing for *Orthoebolavirus* infection, by experienced personnel who have received specific training on the use of the OraQuick™ Ebola 2.0 Rapid Antigen Test, training in the correct use of recommended personal protective equipment (PPE), and expertise in infectious disease diagnostic testing, including the safe handling of clinical specimens potentially containing *Orthoebolavirus*.

PERFORMANCE CHARACTERISTICS

WHOLE BLOOD CLINICAL PERFORMANCE

BUNDIRUGYO POSITIVE CONTRIVED BLOOD SAMPLES

Fifty (50) unique venous whole blood samples including 6 from febrile individuals ($> 100.4^{\circ}\text{F}$) were collected and spiked with inactivated Bundibugyo virus (200706291 Uganda Prototype, lot 815371) to prepare contrived positive samples. Ten (10) samples were spiked at the LoD, fifteen (15) were spiked at $< 2\times$ LoD and twenty-five (25) were spiked in the range of $3\times$ LoD to $5\times$ LoD. Overall positive percent agreement (PPA) for venous whole blood was 100% (95% CI, 93.6%-100%). Results are summarized in the table below.

Sample type	Concentration Tested	Positive percent agreement (n/N)	95% CI
Contrived positive venous whole blood	LoD	100% (10/10)	72.2%-100%
	1.9 x LoD	100% (15/15)	79.6%-100%
	3 x LoD	100% (10/10)	72.2%-100%
	4 x LoD	100% (10/10)	72.2%-100%
	5 x LoD	100% (5/5)	56.6%-100%

n = number in agreement; N = sample size; CI = confidence interval

Six (6) febrile samples were tested as spiked positive samples

SPECIFICITY PERFORMANCE IN NON-ENDEMIC POPULATION

To support venous whole blood and fingerstick whole blood claims, specificity studies were conducted on a total of 245 subjects (166 venous whole blood, 29 febrile venous whole blood and 50 fingerstick whole blood) in the United States. Of the venous whole blood samples, seventy-three (73) unique whole blood samples were also tested for specificity as part of the contrived whole blood samples (prior to spiking). A summary of all results is provided in table below.

Sample type	Negative percent agreement (n/N)	95% CI
Venous whole blood ^a	98.2% (163/166)	94.8%-99.4%
Venous whole blood - febrile ^b	100% (29/29)	89.0%-100%
Fingerstick whole blood	98% (49/50)	88.5%-99.9%

n = number in agreement; N = sample size; CI = confidence interval

^aSixty-six (66) negative samples from the contrived sample study have been incorporated into the total for venous whole blood.

^bSeven (7) febrile samples were tested as part of contrived specimen study.

WHOLE BLOOD ANALYTICAL PERFORMANCE

LIMIT OF DETECTION

A limit of detection (LoD) range-finding study in venous whole blood identified 1.76×10^4 TCID₅₀/mL as the tentative LoD for inactivated Bundibugyo virus (200706291 Uganda Prototype, lot 815371, 3.51×10^6 TCID₅₀/mL). This tentative LoD was confirmed as the LoD by 19/20 replicates testing positive with the same inactivated Bundibugyo virus at this concentration. The range finding study in venous whole blood identified 1.43×10^4 TCID₅₀/mL as the tentative LoD for Ebola Zaire (Mayinga). This tentative LoD was confirmed as the LoD by 19/20 replicated testing positive with the same Mayinga inactivated virus at this concentration. The tentative LoD for Tai Forest (Côte D'Ivoire) was determined to be 1.11×10^3 TCID₅₀/mL in venous whole blood and confirmed in 20/20 replicates. Refer to the table below.

Additionally, a LoD range-finding study was performed using recombinant Mayinga VP40 antigen spiked into pooled venous whole blood. This LoD range-finding study determined the tentative LoD of the OraQuick™ Ebola 2.0 Rapid Antigen Test in venous whole blood with the recombinant antigen to be 0.5 ng/mL. The tentative LoD was confirmed as the LoD by 20/20 replicates testing positive with the same recombinant VP40 antigen at this concentration.

Ebola strain	Concentration tested	Positive
Bundibugyo (200706291 Uganda Prototype)	1.76×10^4 TCID ₅₀ /mL	19/20
Ebola Zaire (Mayinga)	1.43×10^4 TCID ₅₀ /mL	19/20
Tai Forest (Côte D'Ivoire)	1.11×10^3 TCID ₅₀ /mL	20/20

WHOLE BLOOD ANALYTICAL REACTIVITY (INCLUSIVITY)

Analytical reactivity of the OraQuick™ Ebola 2.0 Rapid Antigen Test was evaluated for additional strains of the *Orthoebolavirus*. Testing was performed using negative whole blood as the testing sample matrix. The OraQuick™ Ebola 2.0 Rapid Antigen Test detects all *Orthoebolavirus* strains tested. The results provided below represent the lowest concentration in which all replicates were positive.

Ebola strain	Inactivated or live	Concentration tested	Positive
Sudan (200011676 Gulu)	Inactivated	1.04×10^4 TCID ₅₀ /mL	3/3
	Live	1.04×10^5 TCID ₅₀ /mL	3/3
Reston (119876 Pennsylvania Strain)	Inactivated	3.16×10^7 PFU/mL	5/5

HIGH-DOSE HOOK EFFECT

A high-dose hook effect is a false negative result due to very high concentrations of the target analyte. Antigens at excessive concentrations could quench the specific antibodies in the assay, decreasing complexes between the antibody on the gold conjugate and the biotinylated antibody. Insufficient binding of the gold conjugate could occur at the test line, and the result could be interpreted falsely as a negative result. Testing of live virus stocks at high concentrations of *Orthoebolavirus* species Bundibugyo, Zaire (Mayinga), Sudan and Tai Forest, were evaluated and demonstrated acceptable performance of the OraQuick™ Ebola 2.0 Rapid Antigen Test. A decrease in the visual intensity of the control line was observed; however, the control line remained visible. Refer to table for concentrations tested.

Ebola strain	Concentration tested (TCID ₅₀ /mL)
Bundibugyo (200706291 Uganda Prototype)	3.51 x 10 ⁶
Zaire (Mayinga)	2.85 x 10 ⁷
Sudan (200011676 Gulu)	4.14 x 10 ⁷
Tai Forest (Cote D'Ivoire)	1.00 x 10 ⁶

An analytical study was conducted, and the results demonstrated acceptable performance of the OraQuick™ Ebola 2.0 Rapid Antigen Test when evaluating samples at excessive antigen concentrations. No decrease in visual intensity was observed at the test line for VP40 antigen concentrations up to 50 µg/mL, 100,000 times the established LoD. A decrease in the visual intensity of the control line was observed at higher concentrations; however, the control line remained visible up to 50 µg/mL.

WHOLE BLOOD INTERFERING SUBSTANCES

The OraQuick™ Ebola 2.0 Rapid Antigen Test was evaluated with the following interfering substances present in negative whole blood and whole blood spiked with Zaire Mayinga recombinant antigen (rAg) at ≤ 3x LoD, in order to assess their potential effect on the assay performance. If the substance interfered with the performance of the OraQuick™ Ebola 2.0 Rapid Antigen Test, additional dilutions were prepared to determine at what concentration the substance did not cause interference.

Interfering substance	Target testing concentration	N	Negative samples		Positive rAg-spiked samples	
			Positive	Negative	Positive	Negative
Unconjugated bilirubin	25 mg/dL	3	0/3	3/3	3/3	0/3
Serum protein	5 g/dL	3	0/3	3/3	3/3	0/3
Rheumatoid factor	544 IU/mL	3	2/3	1/3	3/3	0/3
Rheumatoid factor – dilution 1	272 IU/mL	3	0/3	3/3	Not tested	
Cholesterol	13 mmol/L	3	0/3	3/3	3/3	0/3
HAMA	1016 ng/mL	3	0/3	3/3	3/3	0/3
Hemoglobin (hemolyzed blood)	20 g/dL	3	0/3	3/3	1/3	2/3
Hemoglobin (hemolyzed blood) - dilution 1	10 g/dL	3	Not tested		3/3	0/3
ANA	≥ 8.0 AI	3	0/3	3/3	3/3	0/3
Chloroquine	18 µg/mL	3	0/3	3/3	3/3	0/3
Proguanil/Atovaquone	2.8/89.1 µg/mL	3	0/3	3/3	3/3	0/3
Ibuprofen	2425 µmol/L	3	0/3	3/3	3/3	0/3
Acetaminophen	1324 µmol/L	3	0/3	3/3	3/3	0/3
Quinine	148 µmol/L	3	0/3	3/3	3/3	0/3
Amoxicillin	206 µmol/L	3	0/3	3/3	3/3	0/3
Tetracycline	34 µmol/L	3	0/3	3/3	3/3	0/3
Ribavirin	11.04 µg/mL	3	0/3	3/3	3/3	0/3
Biotin	3.6 µg/mL	3	0/5	5/5	5/5	0/5

There was no interference from the exogenous substances tested, including pain relievers, antivirals, antimalarials, or antibiotics. From the endogenous substances tested, false positives were seen with rheumatoid factor (RF) samples at 544 IU/mL; however, rheumatoid factor (RF) samples at 272 IU/mL did not cause interference. False negatives were seen with hemoglobin (hemolyzed blood) samples at 20 g/dL; however, hemoglobin (hemolyzed blood) samples at 10 g/dL did not cause interference. The remaining endogenous substances including bilirubin, cholesterol, serum protein, HAMA (human anti-mouse antibodies), and ANA (antinuclear antibody) did not cause interference with the performance of the OraQuick™ Ebola 2.0 Rapid Antigen Test.

WHOLE BLOOD CROSS-REACTIVITY

The cross-reactivity of the OraQuick™ Ebola 2.0 Rapid Antigen Test evaluated the impact of various virus, bacteria, and parasitic pathogens. In this study, multiple replicates were tested with the pathogens spiked into whole blood at the concentrations listed in table. None of the tested organisms produced false positive results in the OraQuick™ Ebola 2.0 Rapid Antigen Test at the tested concentrations.

Organism	Tested concentration	Units	Positive	Negative
<i>Plasmodium falciparum</i>	unknown	N/A	0/5	5/5
<i>Plasmodium vivax</i>	unknown	N/A	0/5	5/5
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Typhi	2.0E+06	CFU/mL	0/5	5/5
<i>Salmonella enterica</i> subsp. <i>enterica</i> , serovar Choleraesuis	8.4E+07	CFU/mL	0/5	5/5
Dengue virus serotype 1	5.00E+06	TCID ₅₀ /mL	0/5	5/5
Dengue virus serotype 2	8.9E+04	TCID ₅₀ /mL	0/5	5/5
Dengue virus serotype 3	1.6E+06	TCID ₅₀ /mL	0/5	5/5
Dengue virus serotype 4	2.8E+07	TCID ₅₀ /mL	0/5	5/5
Measles virus	8.9E+04	TCID ₅₀ /mL	0/5	5/5
<i>Leptospira biflexa</i> , serotype Patoc	unknown	N/A	0/5	5/5
Simian rotavirus A	1.6E+07	TCID ₅₀ /mL	0/5	5/5
Adenovirus type 4	2.8E+06	TCID ₅₀ /mL	0/5	5/5
<i>Escherichia coli</i> serotype O157:H7	2.30E+07	CFU/mL	0/5	5/5
<i>Shigella dysenteriae</i>	>2E+03	CFU/mL	0/5	5/5
<i>Clostridioides difficile</i>	2.8E+08	CFU/mL	0/5	5/5
<i>Rickettsia africae</i>	1.6E+05	TCID ₅₀ /mL	0/5	5/5
<i>Rickettsia australis</i>	8.9E+05	TCID ₅₀ /mL	0/5	5/5
SARS-CoV-2	1.4E+05	TCID ₅₀ /mL	0/5	5/5
Marburg Ravn	3.39E+05	TCID ₅₀ /mL	0/5	5/5
Marburg Musoke	3.16E+06	PFU/mL	0/5	5/5
Marburg Voegelé *	3.6E+06	PFU/mL	0/5	5/5
Marburg Angola	Unknown	N/A	0/5	5/5
Crimean-Congo hemorrhagic fever virus	5.6E+04	TCID ₅₀ /mL	0/5	5/5
Lassa Josiah	3.16E+05	PFU/mL	0/5	5/5
Rift Valley fever virus	5.43E+06	PFU/mL	0/5	5/5
Yellow fever virus, vaccine strain 17D	1.6E+06	TCID ₅₀ /mL	0/5	5/5
Influenza virus A	3.0E+06	CEID ₅₀ /mL	0/5	5/5
Influenza virus B	2.8E+06	TCID ₅₀ /mL	0/5	5/5
RSV	1.6E+05	TCID ₅₀ /mL	0/5	5/5
Chikungunya	8.9E+06	TCID ₅₀ /mL	0/15	15/15
<i>Coxiella burnetii</i>	1.6E+04	TCID ₅₀ /mL	0/15	15/15
Enterovirus	8.9E+05	TCID ₅₀ /mL	0/5	5/5
<i>Borrelia recurrentis</i>	unknown	N/A	0/5	5/5
<i>Pseudomonas aeruginosa</i>	1.12E+08	CFU/mL	0/5	5/5
<i>Vibrio cholerae</i>	>2E+03	CFU/mL	0/5	5/5
<i>Streptococcus pneumoniae</i>	>1E+03	CFU/mL	0/5	5/5
<i>Hemophilus influenzae</i>	1.8E+07	CFU/mL	0/5	5/5
<i>Neisseria meningitidis</i>	7.2E+05	CFU/mL	0/5	5/5
<i>Yersinia pestis</i>	2.9E+06	CFU/mL	0/5	5/5
<i>Yersinia enterocolitica</i>	1.4E+08	CFU/mL	0/5	5/5
<i>Rickettsia prowazekii</i>	1.6E+03	TCID ₅₀ /mL	0/15	15/15
<i>Rickettsia conorii</i>	2.8E+06	TCID ₅₀ /mL	0/15	15/15
<i>Trypanosoma cruzi</i>	unknown	N/A	0/5	5/5

* Marburg Voegelé is an isolate from a patient infected during the 1967 Marburg outbreak.

CADAVERIC ORAL FLUID

CADAVERIC ORAL FLUID CLINICAL PERFORMANCE

BUNDIBUGYO POSITIVE CONTRIVED ORAL FLUID SAMPLES

Fifty (50) unique oral fluid samples from live donors were collected and spiked with Bundibugyo inactivated virus (200706291 Uganda Prototype, lot 815371) to prepare contrived positive samples as a surrogate for direct collect cadaveric oral fluid samples. Ten (10) samples were spiked at the LoD, fifteen (15) were spiked at < 2x LoD and twenty (25) were spiked in the range of 3x LoD to 5x LoD. Overall positive percent agreement (PPA) for oral fluid from live donors was 100% (95% CI, 93.6%-100%). Results are summarized in table.

Sample type	Concentration Tested	Positive percent agreement (n/N)	95% CI
Contrived positive oral fluid	LoD	100% (10/10)	72.2%-100%
	1.9 x LoD	100% (15/15)	79.6%-100%
	3 x LoD	100% (10/10)	72.2%-100%
	4 x LoD	100% (10/10)	72.2%-100%
	5 x LoD	100% (5/5)	56.6%-100%

n = number in agreement; N = sample size; CI = confidence interval

SPECIFICITY PERFORMANCE IN NON-ENDEMIC POPULATION

To support cadaveric oral fluid claims, specificity studies were conducted on 108 living donors (oral fluid) and 80 cadavers in the United States. Oral fluid from living donors was tested as a surrogate for cadaveric oral fluid. For living donors, the specificity demonstrated was 100% (108/108; 95% CI, 96.6%-100%) for direct collection. For cadaver samples, the direct collection method demonstrated 100% (80/80; 95% CI, 95.5%-100%). A summary of all results is provided in table.

Sample type	Negative percent agreement (n/N)	95% CI
Oral fluid — live donors ^a	100% (108/108)	96.6%-100%
Cadaveric oral fluid	100% (80/80)	95.5%-100%

n = number in agreement; N = sample size; CI = confidence interval

^aFifty-eight (58) negative samples were collected and tested as part of the contrived sample study.

CADAVERIC ORAL FLUID ANALYTICAL PERFORMANCE

LIMIT OF DETECTION

A limit of detection (LoD) range-finding study was conducted in oral fluid collected from living donors spiked onto the collection pad of the OraQuick™ Ebola 2.0 Rapid Antigen Test. The range-finding study identified 8.78×10^3 TCID₅₀/mL as the tentative LoD for Bundibugyo inactivated virus (200706291 Uganda Prototype, lot 815371, 3.51×10^6 TCID₅₀/mL). This tentative LoD was confirmed as the LoD by 19/20 replicates testing positive with the same Bundibugyo inactivated virus at this concentration.

CROSS-REACTIVITY — CADAVERIC ORAL FLUID

The cross-reactivity of the OraQuick™ Ebola 2.0 Rapid Antigen Test was evaluated by testing additional viral and bacterial pathogens. In this study, five (5) replicates were tested with the pathogens spiked into Ebola negative oral fluid (saliva as a surrogate) at the concentrations listed in table. None of the tested organisms produced false positive results in the OraQuick™ Ebola 2.0 Rapid Antigen Test.

Organism	Tested concentration	Units	Positive	Negative
<i>Moraxella (Branhamella) catarrhalis</i>	3.8E+06	CFU/mL	0/5	5/5
<i>Corynebacterium diphtheriae</i>	3.2E+07	CFU/mL	0/5	5/5
<i>Nocardia spp.</i>	Unknown	N/A	0/5	5/5
<i>Hoylella (Bacteroides) oralis</i>	>2E+03	CFU/mL	0/5	5/5
<i>Chlamydomydia pneumoniae</i>	2.9E+07	IFU/mL	0/5	5/5
<i>Mycoplasma pneumoniae</i>	2.6E+07	CFU/mL	0/5	5/5
<i>Bordetella pertussis</i>	8.0E+07	CFU/mL	0/5	5/5
<i>Clostridioides difficile</i>	2.8E+08	CFU/mL	0/5	5/5
Cytomegalovirus (human herpesvirus 5)	2.8E+04	TCID ₅₀ /mL	0/5	5/5
Epstein-Barr virus (human gammaherpes virus 4)	Unknown	N/A	0/5	5/5

Additional cross-reactivity was performed for the OraQuick™ Ebola 2.0 Rapid Antigen Test for use with whole blood and no cross-reactivity was observed at the tested concentrations. Refer to the Whole Blood Cross-Reactivity testing section for the list of organisms.

BIBLIOGRAPHY

1. CDC, Universal Precautions for Prevention of Transmission of Human Immunodeficiency Virus, Hepatitis B Virus, and other Bloodborne Pathogens in Health-Care Settings. MMWR 1988; 37(24):377-388.
2. Infection Control for Viral Haemorrhagic Fevers in the African Health Care Setting. Centers for Disease Control and Prevention and World Health Organization. Atlanta, Centers for Disease Control and Prevention, 1998: 1-198. WHO, September 2014: Ref: WHO/HIS/SDS/2014.4
3. *EMERGENCY GUIDELINE Implementation and management of contact tracing for Ebola virus disease* (<https://iris.who.int/server/api/core/bitstreams/1be2b0a4-ecaa-442d-af2b-709dae0dae58/content>) issued by the World Health Organization (WHO) and the Centers for Disease Control (CDC).

EXPLANATION OF SYMBOLS

 LOT	Batch Code	 IVD	In Vitro Diagnostic Medical Device
 REF	Catalog Number		Manufacturer
	Caution, Consult Accompanying Documents	 PN	Part Number
	Use By		Temperature Limitation
		 R_x	Medical Prescription

For Technical or Customer Service within the United States, 1-800-ORASURE (1-800-672-7873).
For customers outside the United States, phone +(001) 610-882-1820 or go to www.orasure.com.

220 East First Street
Bethlehem, PA 18015 United States
ORASURE (1-800-672-7873) • (610) 882-1820
www.orasure.com

© 2026 OraSure Technologies, Inc.
This product is covered by one or more patents.
OraQuick is a trademark of OraSure Technologies, Inc., registered in various jurisdictions.

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BIBLIOGRAPHIE

1. CDC, Universal Precautions for Prevention of Transmission of Human Immunodeficiency Virus, Hepatitis B Virus, and other Bloodborne Pathogens in Health-Care Settings. MMWR 1988; 37(24):377-388.
2. Infection Control for Viral Haemorrhagic Fevers in the African Health Care Setting. Centers for Disease Control and Prevention et Organisation mondiale de la Santé. Atlanta, Centers for Disease Control and Prevention, 1998 : 1-198. OMS, Septembre 2014 : Réf : WHO/HIS/SDS/2014.4
3. *DIRECTIVE RELATIVES AUX SITUATIONS D'URGENCE Mise en œuvre et gestion de la recherche de contacts pour la maladie à virus Ebola* (<https://iris.who.int/server/api/core/bitstreams/1be2b0a4-ecaa-442d-af2b-709dae0dae58/content>) publiée par l'Organisation mondiale de la Santé (OMS) et les Centers for Disease Control (CDC).

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OraQuick est une marque d'OraSure Technologies, Inc., déposée dans différentes juridictions.