For use under Emergency Use Authorization (EUA) only.
For in vitro diagnostic use.

The OraQuick™ Ebola 2.0 Rapid Antigen Test is an in vitro diagnostic single-use immunoassay for the qualitative detection of antigens from viruses within the *Orthoebolavirus* genus but does not differentiate between these viruses.

This product has not been U.S. FDA cleared or approved, but has been authorized for emergency use by U.S. FDA under an EUA for use by the authorized laboratory.

This product has been authorized only for the detection of antigens from orthoebolaviruses including *Orthoebolavirus zairense*, *Orthoebolavirus bundibugyoense*, *Orthoebolavirus sudanense*, and *Orthoebolavirus taiense*, not for any other viruses or pathogens.

The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection of Ebola virus under section 564(b)(1) of the Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.

This Quick Reference Guide is designed to provide a basic guideline for testing. Please refer to the package insert for complete information on proper collection of samples and procedural use of the OraQuick™ Ebola 2.0 Rapid Antigen Test. Deviation from these procedures may result in inaccurate results.

POUCH CONTENTS: 1 developer vial, 1 device, and 1 absorbent packet

STEP 1. – Specimen Collection ^{1,2,3}

Using the flat pad of the test device, swab the upper and lower gum line as shown in the photos below. Swab the cadaver using the flat pad of the device. After collection, go to Step 2.

Collect sample from gum line.
Swab lower and upper gum once.

OR
Collect sample from one inner
cheek. Swab once.

OR
Collect sample from the soft palate.
Swab once.

STEP 2. – Perform the Test

Insert the collection pad of the device all the way into the developer vial.
Make sure that the collection pad touches the bottom of the vial.
The result window on the device should be facing you.

Set a timer for 30 minutes
and start timing the test.

DO NOT remove the device from the
developer vial while the test is running.

Fluid will appear and travel up the result window.
The fluid may gradually lighten as the test develops.

STEP 3. – Read Results

Read the results in a fully lit area. Read results between 30 and 40 minutes to obtain an accurate result. DO NOT read before 30 minutes or after 40 minutes.

To report a device problem, contact OraSure Technologies Customer Service (1-800-orasure) within the U.S. or +(001) 610-882-1820 for outside the U.S. You can also go to www.orasure.com/inquiry.html.

Negative: Line in C zone

Line in C zone

- A negative test result means that Ebola antigens were not detected in the specimen.
- The test result is interpreted as Ebola antigens were not detected. The cadaver is negative for Ebola antigens.

Invalid: Repeat with a fresh sample

No lines
present

Purple background
obscures results

Any partial line
at C or T zone

Positive: Lines in T zones

Example of positive results.

Line in T zone

Line in T zone

Faint line in T zone

No line in C zone

- A positive test result means that Ebola antigens have been detected in the specimen. The cadaver is presumed positive for Ebola.
- Individuals with a positive result with the OraQuick™ Ebola 2.0 Rapid Antigen Test should, in accordance with CDC and WHO recommendations, be subjected to a safe and dignified burial procedure. Contacts of an Ebola-positive cadaver should be identified and monitored.
- In some cases, the purple line in the C Zone may not be present or may be very faint if there are high levels of virus in the sample.

Order Information and Description

1001-0708 25-ct kit
1001-0709 kit controls

- See "Universal Precautions," CDC, MMWR, 1988; 37(24):377-388.
- Centers for Disease Control and Prevention and World Health Organization. Infection Control for Viral Haemorrhagic Fevers in the African Health Care Setting. Atlanta, Centers for Disease Control and Prevention, 1998:1-198.
- Information for Health Care Workers in the United States: <http://www.cdc.gov/vhf/ebola/healthcare-us/index.html>

© 2026 OraSure Technologies, Inc. OraQuick Ebola is a trademark of OraSure Technologies, Inc., registered in various jurisdictions.
3001-3962_EUA v2/07-26

**OraSure
Technologies, Inc.**
220 East First Street
Bethlehem, PA 18015 USA
phone: 800.ORASURE
web: www.orasure.com
Made in the USA

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