For use under Emergency Use Authorization (EUA) only.
For in vitro diagnostic use.

The OraQuick™ Ebola 2.0 Rapid Antigen Test is an in vitro diagnostic single-use immunoassay for the qualitative detection of antigens from viruses within the *Orthoebolavirus* genus but does not differentiate between these viruses.

This product has not been U.S. FDA cleared or approved, but has been authorized for emergency use by U.S. FDA under an EUA for use by the authorized laboratory.

This product has been authorized only for the detection of antigens from orthoebolaviruses including *Orthoebolavirus zairense*, *Orthoebolavirus bundibugyoense*, *Orthoebolavirus sudanense*, and *Orthoebolavirus taiense*, not for any other viruses or pathogens.

The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection of Ebola virus under section 564(b)(1) of the Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.

This Quick Reference Guide is designed to provide a basic guideline for testing. Please refer to the package insert for complete information on proper collection of samples and procedural use of the OraQuick™ Ebola 2.0 Rapid Antigen Test. Deviation from these procedures may result in inaccurate results.

POUCH CONTENTS: 1 developer vial, 1 device, and 1 absorbent packet

STEP 1a. – Collect the Specimen^{1,2,3} Fingerstick Whole Blood

Cleanse finger. Air dry.
Puncture with lancet.

Capillary blood tube

Handle

Wipe away first drop of blood. Fill capillary blood tube to black lines.

NOTE: Touch the tip of the capillary blood tube to the drop of blood. Filling is automatic; never squeeze the capillary blood tube while sampling. To prevent bubbles, the tip of the capillary blood tube must maintain uninterrupted contact with the blood sample until filled.

If using a capillary blood tube, make sure it is filled with blood up to the indicator line and no bubbles are present.

STEP 1b. – Collect the Specimen^{1,2,3} Venipuncture Whole Blood

Collect blood using standard phlebotomy procedures.

- Prior to testing, mix blood tube gently by inversion several times to ensure homogenous specimen.
- Remove tube cap.
- Whole blood may be stored for 6 days refrigerated at 2°C-8°C (36°F-46°F).

Using a calibrated pipette, pipette 50 µL of sample and inspect the tip for bubbles. If a bubble is present, discard the collected sample and obtain a new sample.

OR

STEP 2. – Perform the Test Fingerstick Whole Blood and Venipuncture Whole Blood

Deposit the blood sample by compressing the bulb on the capillary blood tube, or from the calibrated pipette onto the collection pad at the location shown in the photo. Samples MUST only be applied directly to the device and should not be pipetted into the developer vial. **Note: Pipetting the sample at a location other than that shown in the photo could result in failure to detect a low positive sample.**

Insert the collection pad of the device all the way into the developer vial. Make sure that the collection pad touches the bottom of the vial. The result window on the device should be facing you.

Set a timer for 30 minutes and start timing the test.

Fluid will appear and travel up the result window. The fluid may gradually lighten as the test develops.

DO NOT remove the device from the developer vial while the test is running.

STEP 3. – Read Results

Read the results in a fully lit area. Read results between 30 and 40 minutes to obtain an accurate result. DO NOT read before 30 minutes or after 40 minutes.

To report a device problem, contact OraSure Technologies Customer Service (1-800-orasure) within the U.S. or +(001) 610-882-1820 for outside the U.S. You can also go to www.orasure.com/inquiry.html.

Negative: Line in C zone

Line in C zone.

- A negative test result means that Ebola antigens were not detected in the specimen.
- The test result is interpreted as Ebola antigens were not detected. The patient is presumed negative for Ebola antigens.

Invalid: Repeat test

No lines present

Background obscures results

Any partial line at C or T zone

Positive: Lines in T zones

Examples of positive results.

Line in T zone

Line in T zone

Faint line in T zone

No line in C zone

- A positive test result means that Ebola antigens have been detected in the specimen. The patient is presumed positive for Ebola antigens. Follow up testing with PCR to confirm test result.
- Some pink may remain visible in the background as a result of the whole blood specimen.
- Immediately follow the CDC guidelines for at-risk exposure patient management.
- In some cases, the purple line in the C zone may not be present or may be very faint if there are high levels of virus in the sample.

Order Information and Description

1001-0708 25 Count kit
1001-0709 Kit Controls

1 See "Universal Precautions," CDC, MMWR, 1988; 37(24):377-388.

2 Centers for Disease Control and Prevention and World Health Organization, Infection Control for Viral Haemorrhagic Fevers in the African Health Care Setting. Atlanta, Centers for Disease Control and Prevention, 1998:1-198.

3 Information for Health Care Workers in the United States: <https://www.cdc.gov/ebola/about/index.html>

© 2026 OraSure Technologies, Inc. OraQuick Ebola is a trademark of OraSure Technologies, Inc., registered in various jurisdictions.
3001-3963_EUA v2/07-26

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