



AMABWIRIZA Y'IMIKORESHEREZE

Ugomba gukurikiza neza amabwiriza yo kwisuzuma virusi itera SIDA kugirango ubone igisubizo nyacyo. Ntabwo ugomba kugira icyo urya cyangwa umywa nibura mu minota 15 ibanziriza kwipima. Kandi ntabwo ugomba kwoza mu kanwa ukoresheje umiti y' amenyi mu minota 30 ibanziriza kwipima.

IKITONDERWA: Niba ufata imiti igabanya ubukana bwa VIRUSI itera SIDA ushobora kubona igisubizo kitari cyo.



REBA AMABWIRIZA

www.oraquickhivselftest.com

MU KINYARWANDA

ORAQUICK®

UBURYO BWIHUSE BWO KWISUZUMA VIRUSI ITERA SIDA BWA

UKO BAKOresha UBURYO BWIHUSE BWO KWISUZUMA VIRUSI ITERA SIDA UKOreshEJE ORAQUICK®



URAKENERA UBURYO BWO KUBARA IGIHE IKI KIZAMINI KIMARA



Muri aka gashashi harimo: **igikoresho cyo kwisuzumisha, icyo bagiterakaho bapima hamwe n'amabwiriza ku mikoreshereze yacyo.**



Aka gashashi karimo utwumba tubiri.



Ca agashashi mu kumba karimo agacupa.



Gafungure.




BIRABUJIJWE kumena amazi arimo. **NTANYWEBWA.**



Teraka agacupa mu mwanya gaterekwamo.



Ca agashashi uvanemo igikoresho cyo kwipimisha kirimo. **NTUFATE** agatwe k'igikoresho cyo kwipimisha n' intoki zawe. **Nturye** cyangwa ngo umire umuti.



Agatwe k'igikoresho cyo kwipimisha gakoze ukomeje kw'ishinya yo hejuru usa nk'uyihanagura (reba ifoto ya mbere) kandi ubikore no kw'ishinya yo hasi (reba ifoto ya kabiri).



Shyira agatwe k'igikoresho cyo kwipimisha mu gacupa karimo wa muti ugeze neza ku ndiba.




REKA KAMAREMO IMINOTA MAKUMYABIRI (20) mbere yo kukavanamo ngo usome ibisubizo. NTUGOMBA gusoma ibisubizo nyuma y'iminota 40.

UKO BASOMA IBISUBIZO



Somera ibisubizo ahantu hacanye neza

IGISUBIZO KEREKANA KO USHOBORA KUBA YARANDUYE VIRUSI ITERA SIDA



Nk'ikimenyetso cy'uko umuntu yanduye VIRUSI itera SIDA, hagaragara imirongo ibiri yuzuye, nubwo yaba itagaragara neza. Igihe ubonye iki kimenyetso, ihutire kujya kongera kwisuzumishiriza mu kigonderabuzima.



Uhutire...

Njya ku kigonderabuzima kikwegereye gisuzuma VIRUSI itera SIDA

IGISUBIZO KITARI CYO



Iyo nta murongo uri kuri ca "C" (n'iyu haba hari umurongo kuri ta "T") cyangwa ukabona akadirishya gasa nk'agaturu gatuma biragora kubona igisubizo, bigaragaza ko igisubizo kitari cyo.

Biragusaba kongera kwisuzuma ukoresheje ikindi gikoresho cyo kwipima cya OraQuick.



Ibizamini ntibyafashwe neza.

Njya ku kigonderabuzima kikwegere gisuzuma VIRUSI itera SIDA wongere wipimisha.

IGISUBIZO KEREKANA KO USHOBORA KUBA UTARANDUYE VIRUSI ITERA SIDA

NIBA USOMYE IBISUBIZO MBERE Y'IMINOTA 20, USHOBORA KUBONA IGISUBIZO KITARI CYO



Iyo hari UMURONGO UMWE KURI CA "C" ariko ntubone uri kuri TA "T", bivuga ko UTANDUYE VIRUSI itera SIDA.

Isuzumisha ahantu abantu basanzwe bisuzumishiriza. Niba hari aho ushobora kuba ukeka kuba waranduriye Virus ya HIV itera SIDA, ongera wisuzumisha nyuma y'amezi 3.

IBISUBIZO NTIBYZEWE

Ntuzi cyangwa ntizizwe igisubizo cyawe.

Njya ku kigonderabuzima kikwegere gisuzuma VIRUSI itera SIDA wongere wipimisha.

BIJUGUNYE

Vanamo akambi wakoresha wisuzuma, upfundikire agacupa hanyuma ujugunye ibirimo byose ahabugenewe.

